
ENFORCE ACTIVE FOOT

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www.enforcemed.pl

Contents

Product description.....	8
Construction	8
Technical specifications	8
Application	9
Prosthetic foot selection	9
Contraindications.....	9
Conditions of use and storage	10
Care and maintenance of the product	10
Safety instructions	10
Responsibility	11
Warranty	11
Declaration of conformity	11
Disposal	11

- Date of the last update of the document: 08.10.2024 (IAF.10.2024)
- Before using the product, read this document and follow the safety instructions.
- Any incidents with negative health effects related to the product should be immediately reported to the manufacturer by sending a report via e-mail to: complaints@enforcemed.pl. Include such information as: name, date of occurrence and description of the situation, serial number of the product found on the label.
- Inform medical personnel of any changes in weight that make it necessary to change the foot (with reference to the table under prosthetic foot selection) or activity level.

Product description


The term device refers to the composite prosthetic foot - ENFORCE Active Foot manufactured by ENFORCE Medical Technologies LLC.

The ENFORCE Active Foot prosthesis consists of two composite components along with a male pyramid adapter. The prosthesis is designed for lower limb amputees with activity level II-IV. It is adapted to move on different types of ground along with variable speed.

Construction

The product consists of the following components:

1. Composite upper component.
2. Composite upper component.
3. Flexible heel wedge.
4. Foot adapter.

Technical specifications

Component weight (23 cm, P5)	400 g
Maximum patient weight	130 kg
Activity level	2 - 4
Available sizes	37 - 45 EU

Application

ENFORCE Active Foot is a non-invasive prosthetic device that replaces the function of the foot, designed for multiple use by a single patient. It has applications for both low and high activity levels. The prosthetic foot, thanks to the Forefoot Dynamic Properties FDP™ system, provides high energy return, a slit in the toe area for greater stability and good ground grip during gait.

This product is compatible with the ENFORCE modular lower limb prosthesis system, functionality with components from other manufacturers has not been tested.

Prosthetic foot selection

Caution::

The prosthetic foot is selected individually based on the patient's current weight, if there are any changes, the medical staff should be informed immediately.

The ENFORCE Active Foot is intended for patients with a mobility level of 2,3 or 4 (for patients who can move freely at home and outdoors). Use beyond our recommendations should be subject to a detailed and individual evaluation by the medical personnel responsible for the prosthesis and familiar with the patient's history.

The device can be configured and fitted to the patient only by authorized persons (trained medical personnel).

The following table provides information on the proper selection of the prosthetic foot for the patient's weight.

Body weight (kg)	<50	<90	<130
Type	P3	P4	P5

Contraindications

The device is not intended for patients with activity level 1 or with a body weight exceeding 130 kg.

Conditions of use and storage

The product is weatherproof, meaning that it can be used in wet or humid environments and is resistant to splashing with fresh water (such as rain). The product should be wiped dry with a clean cloth after contact with fresh water or moisture.

In case of contact with chlorinated water, chemicals, sand, dust or dirt, clean the product immediately with fresh water and dry.

Make sure the product is completely dry before reusing it.

The product's usable temperature range: - 50°C to +100° C

Storage and transportation: temperature range -50°C to +100°C, relative humidity up to 90%.

Care and maintenance of the product

The product can be cleaned with clean fresh water mixed with pH-neutral soap, finally wiped with a dry cloth and allowed to dry completely. Do not use caustic cleaners.

After the first 60 days of active use, the product should be visually inspected. In addition, it is recommended to conduct annual safety inspections performed only by trained personnel.

Safety instructions

Caution::

When the product shows signs of damage, stop using it immediately.

Signs of change or loss of utility function:

- Changing gait characteristics.
- Changing the position of components.
- Appearing of abnormal noises in prosthesis components.

Important safety information which must be followed carefully:

- Use the product only during the approved service life designated by the manufacturer (3 years).
- Submit the product to an annual safety inspection.
- If damage to the product is suspected, check its operation and serviceability.
- If the product is damaged, stop using it immediately.

- Use only under permitted conditions and avoid corrosive environments.
- The product should be inspected for damage if it has been used in unapproved environmental conditions.
- Clean the contaminated product according to the instructions.
- When descending stairs and steep surfaces, use a handrail if available.
- The device is not suitable for extreme sports, running or cycle racing. Only occasionally bike rides are allowed. Any such activities undertaken are done so completely at the users' own risk.
- Do not place body parts between the composite elements of the product, as this may result in loss of skin continuity.
- The product should be used with a dedicated cosmetic cover.
- Be especially careful when walking on uneven or slippery surfaces.

Responsibility

The manufacturer (ENFORCE Medical Technologies LLC) is liable if the product is used in accordance with the recommendations contained in this document.

Warranty

ENFORCE Medical Technologies provides a 24-month warranty on the product. Please refer to the product warranty card for details. Changes or modifications not expressly approved may void the warranty.

Declaration of conformity

The product meets the requirements of Regulation (EU) 2017/745 of the European Parliament and of the Council of April 5, 2017 on medical devices.

Disposal

Dispose of the product and packaging in accordance with applicable local or national environmental regulations.

Medical Device

The product meets the requirements of Regulation (EU) 2017/745 of the European Parliament and of the Council of April 5, 2017 on medical devices

ENFORCE Medical Technologies Sp. z o.o.,
ul. Szczęśliwowska 45A 61-626 Poznań, Poland.
tel.: 61 111 00 69,
email: biuro@enforcemed.pl, www.enforcemed.pl

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