# ENFORCE BIONIC FOOT

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**Table of contents**

1. Product description .....	38
2. Operation instructions .....	40
3. Mobile application.....	49
4. Conditions of use and storage.....	66
5. Maintenance and care .....	67
6. Legal and safety notices .....	67
7. Responsibility .....	68
8. Warranty information .....	68
10. Wireless modules.....	70
11. Declaration of conformity (CE).....	70
12. Disposal .....	70

- Date of last update of the document: 19.03.2024 (IBF.03.2024)
- Before using the product, read this document and observe the safety instructions.
- Any incidents with a negative impact on health related to the product should be reported immediately to the manufacturer by sending an email to [complaints@enforcemed.pl](mailto:complaints@enforcemed.pl). Please include information such as your name, the date of the occurrence, a description of the situation, and the serial number of the product found on the label on the packaging. Additionally, the incident should be reported to the relevant authority responsible for the place of residence.
- Any changes in weight exceeding four kilograms should be communicated to medical staff.

## 1. Product description

The term device refers to the electronic ankle prosthesis - ENFORCE Bionic Foot manufactured by ENFORCE Medical Technologies LLC. The product assists walking with active damping of foot movements in the sagittal plane, controlled by a microprocessor. The device takes real-time measurements from built-in sensors to change the state of the magneto-rheological damper according to the current gait phase. The device takes measurements at a frequency of 100Hz, which ensures low latency of the control algorithm - responding dynamically to changes in the state of the prosthesis. The damping characteristics are dependent on the user's gait pattern - adapting to the range of changes in sole and dorsiflexion. It is possible to adjust the prosthesis to the individual needs of the user by trained medical personnel.

The ENFORCE Prosthetic app is used for operating the device, providing the following functionalities:

- Connecting to the device via Bluetooth.
- Reading the battery charge status.
- Reading and adjusting the damper power.
- Changing the operating mode.
- Collecting statistics of steps taken on the enabled prosthesis.
- Service notifications.
- Software updates.

The application is designed to adjust the parameters of the electronic prosthesis according to the current footwear, compensating for the influence of its shape on the performance of the control algorithm.

The device is equipped with LEDs arranged in a ring, serving to convey visual messages to the user regarding the prosthesis's status without the need for connecting to an app. The LED indicator provides information on:

- Percentage of battery charge.
- Current battery level during charging.
- Full charge status when connected for charging.
- Connection and disconnection of the device with the app via Bluetooth.
- On and off functions.
- Progress of a software update.
- Errors in the operation of the prosthesis.

## 1.1 Product specifications

1. Range sensor
2. Ring of LEDs
3. Charging connector
4. Magnetorheological damping unit
5. Wiring harness


Figure 1: Construction of the product.

## 1.2 Product applications

The ENFORCE Bionic Foot is a non-invasive prosthetic device designed to replace foot function and is suitable for multiple uses by a single patient. It is applicable for both Activity Level 2 and Activity Level 3, addressing the needs of patients with limited or no mobility, especially in open areas.

This product seamlessly integrates with the ENFORCE modular lower limb prosthesis system. For further details regarding compatibility with products from other manufacturers, please contact the ENFORCE staff.

### 1.3 Prosthetic foot selection

**Attention:**

The prosthetic foot is individually selected based on the patient's current weight. If there is a change in weight, the patient should promptly notify medical personnel to ensure proper adjustment and continued optimal use of the prosthetic

**Attention:**

The device can only be configured and fitted to the patient by authorized individuals, specifically trained medical personnel.

Any use of the device beyond the recommendations provided in this manual should undergo a thorough and individual evaluation by the medical personnel responsible for the prosthesis, who should be familiar with the patient's medical history.

The table below contains information on the appropriate selection of a prosthetic foot for the patient's weight.

Body weight (kg)	<60	<80	<100
Type	P3	P4	P5

### 1.4 Limitations

The product is not intended for use in active sports.

## 2. Operation instructions

### 2.1 ON/OFF

The device can be activated when worn on a limb by taking a step with a firm heel strike against the ground. Alternatively, it can be activated while sitting by lightly striking the ground from a height of about 5 cm. Activation is indicated by an LED ring, and an animation is displayed as shown in the *Table 1*.

The device automatically turns off after 2 minutes of inactivity. Shutdown is indicated by the LED ring, and an animation is displayed as shown in the *Table 1*.

## 2.2 Battery

### 2.2.1 Charging

**Attention:**

Before starting the charging process, it is necessary to remove the product and place it on a stable surface in an upright position.

**Attention:**

The device and its charger must not be used by individuals with pacemakers or implantable cardioverter defibrillators (ICDs).

**Attention:**

The product cannot be charged when the outside temperature exceeds +45°C and at temperatures below 0°C.

### 2.2.2 Charger connection guide

**Attention:**

Before connecting the charger, make sure that there are no foreign objects or moisture in the charging socket and that the charging cable is not damaged.

**Attention:**

Before connecting another charger, ensure its compatibility by checking the parameters against the following list:

- Voltage 5.0-6.2V.
- Minimum current 2.4A.
- USB 2.0 or 3.0

### 2.2.3 Charging process

Review the notes above in sections 2.2.1 to 2.2.3, and then:

- Connect the charging cable to the device.
- Charging has started when the LED Ring displays the current charging status.
- Charging is complete when the entire LED Ring pulses evenly (alternating between on and off every second).

The battery charge level is indicated on the LED Ring. An example animation sequence is shown in the *Table 1*. Charging can be interrupted at any time.

*Example:* Flashing fourth LED indicates charging in the range of 31-40%.

The total charging time is approximately 6 hours with a 12W charger (5V-2.4A). In the case of a deeply discharged battery, a preliminary charging process may begin. This may result in a longer charging process.

#### **2.2.4 Charge level of the prosthesis**

To check the charge level of the prosthesis directly on the LED panel, tilt the upper part of the prosthesis to the maximum (plantar flexion) as shown in *Figure 2*. After 2 seconds, the LED panel will display the current charge level of the prosthesis. Each LED on the panel corresponds to 10% of the battery charge.


Figure 2: Prosthesis position in plantar flexion.

#### **2.2.5 Low battery level**

When the battery level drops below 10%, the device goes into battery saving mode. Therefore, some functions, such as adding new shoes and launching updates, will be automatically disabled, only the option of reducing the damper power and turning off the LED ring remains. This is indicated by the LED ring, one red LED at the top of the LED Ring lights up. The illustration is in the table1.

### **2.3 Calibration**

The calibration process is used to adjust the prosthesis to the footwear. All steps should be performed with the prosthesis and chosen footwear in place, and the description of each calibration step should be followed for its proper execution.

### 2.3.1 Step 1

Stand in a natural position, then press the calibrate button in the mobile application, and remain still until the prosthesis completes the calibration. Moving the prosthesis during this procedure may result in improper functioning.

### 2.3.2 Step 2

Put the tip of the shoe against a wall, for example. Position the ankle joint in neutral, so that the upper part of the prosthesis is perpendicular to the ground. Press the calibrate button, moving the prosthesis during this procedure may cause it to malfunction.

### 2.3.3 Step 3

In this step, perform dorsal flexion of the prosthesis (for example, by performing a squat), the exact positioning of the prosthesis is shown in *Figure 3*. Then, press the calibrate button and, in stillness, wait for the prosthesis to perform calibration. Moving the prosthesis during this procedure may cause it to malfunction.


Figure 3: Prosthesis position in dorsal flexion.

### 2.3.4 Step 4

Stand in a natural position, enter your weight in the app, then press the calibrate button. In stillness, wait for the prosthesis to perform the calibration. Moving the prosthesis during this procedure (changing the load on the prosthesis) may cause it to malfunction.

### 2.3.5 Step 5

In the next stage of calibration, you will need to go through several walking steps. When you are ready, press the START button and walk a minimum of 5 steps. Once done, press the STOP button.

### 2.3.6 Step 6

In the final stage of prosthesis calibration, you will need to perform a squat with the maximum possible depth. After performing the squat, press the calibrate button and, in stillness, wait for the prosthesis to complete the calibration. Moving the prosthesis during this procedure (changing the load on the prosthesis) may cause it to malfunction.

## 2.4 Behavioral modes

### 2.4.1 Walking mode

The correct use of the prosthesis depends on the patient's fitness level. Always adhere to the recommendations provided by the prosthetic staff. The magnetorheological damper supports active toe-off during gait, the damper is switched off when the toes are detached from the ground. In this mode, you can adjust the maximum stiffness of the prosthesis using the mobile app. Refer to section 3.11 for detailed instructions on this process.

**Activation** occurs when every other mode is disabled.

**Deactivation** occurs when one of the available modes is activated via the mobile app.

**Recommendations:** It is advisable to use this mode for the majority of the prosthesis's operating time, as in this mode, the prosthesis automatically adapts to the conditions of gait.

### 2.4.2 Standing assistance mode

This mode is designed to actively support the standing position. When active, the prosthesis actively prevents bending of the ankle joint, irrespective of the ground angle on which the user is standing.

**Activation** automatically during ankle joint immobility sustained for 10 seconds. It is operational exclusively in walking mode and house mode.

**Deactivation** occurs at the onset of walking or with an abrupt rise of the shin.

### 2.4.3 House mode

The mode has been optimized for indoor use and short-distance walking. In this configuration, the prosthesis adjusts the level of support based on gait speed, considering the maximum selected stiffness level in the mobile app.

**Activation/Deactivation** is performed using both the gesture and the application in the Modes panel, as described in the *section 3.10*. To change the state (enable or deactivate the mode) using a gesture (see Figure 4), hold the prosthesis in maximum dorsal flexion for 10 seconds.


Figure 4: Prosthesis position in dorsal flexion.

#### 2.4.4 Bike mode

To facilitate and ensure predictable performance in bike mode, the prosthesis maintains a constant level of damping, matching the resistance of the magnetorheological damper.

**Activation/Deactivation** only through the mobile application in the Modes section described in *point 3.10*.

**Recommendations:** it is recommended to use the mode while cycling.

#### 2.4.5 Standby mode

In this mode, damping is inactive, allowing the foot to fall freely when sitting.

**Activation** when the prosthesis remains immobile for 2 minutes. In this mode, the prosthesis consumes significantly less energy, thereby extending its service life.

**Deactivation** occurs when the foot is lightly stomped or impacted on the ground. The prosthesis is ready for use once the actuation procedure is completed.

### 2.4.6 Sleep mode

The mode allows the prosthesis to be completely switched off. Activating this mode enables the maximum extension of battery life.

**Activation** the mode is possible only through the mobile application in the modes section, according to *point 3.10*.

**Deactivation** the mode is only possible by connecting the charger, before doing so, please refer to the *2.2.1 section* of this manual.

### 2.4.7 Ladder mode

**Attention:**

In this mode, power consumption increases significantly. Use only when necessary.

This mode allows complete stiffening of the prosthesis regardless of external conditions.

**Activation** the mode can only be activated via the mobile application in the modes section, in accordance with the *3.10 point*. After activating the mode in the application, stand in a naturally upright position, as in the first calibration step. Only then will the mode be activated.

**Deactivation** is possible via the mobile application, as well as by setting the prosthesis in a horizontal position, e.g. by bending the knee joint at an angle of 90 degrees.

## 2.5 Product functions overview

### 2.5.1 Bluetooth

Bluetooth functionality is utilized for communication with the mobile app. To extend the prosthesis's working time, we recommend disabling the Bluetooth connection when not actively connecting to the app.

**Activation** is done by placing the prosthesis in a specific position as shown in the illustration 5 (for example, by taking off the prosthesis and rotating the foot upward) for 5 seconds. The activation of Bluetooth is indicated by the steady flashing of half of the LED ring in blue and the other half in green, as shown in the *Table 1*.

**Deactivation** is done similarly. Bluetooth deactivation is signaled by the even pulsating of half of the LED ring in blue and the other half in red, as shown in the *Table 1*.


Figure 5: The position of the prosthesis required to turn on/off Bluetooth.

### 2.5.2 LED indication ring

Action	Animation	Product status
Turn on		After the animation completes, the product is ready to begin walking
Turn off		When the animation completes, the product transitions into standby mode
Connecting to the mobile app		Connecting to the mobile application is underway
Disconnecting with the mobile app		Disconnecting to the mobile application is underway
Connected to the mobile app		The product is connected to a mobile application
Low battery level		Low battery, certain functions disabled. Please charge the product at your earliest convenience
Confirmation		Confirmation of correct change of prosthesis parameters







Action	Animation	Product status
Starting the update		The product update process has started, waiting for the device to join
Update progress		Update in progress, please wait
Battery charging		The product displays the current charging stage
Battery precharging		The product displays the current precharging stage
Battery charge level		The product indicates the current level of charge
Enabling mobile application connectivity		The product signals the activation of Bluetooth connectivity functionality
Disabling mobile application connectivity		The product signals the deactivation of Bluetooth connectivity functionality

Table 1: Summary of LED ring animation.



Action	Animation	Product status
Warning		Warning, contact your prosthetics company
Error		Error, contact the prosthetics company immediately
Critical error		Malfunction, do not use the product. Contact the prosthetics company immediately.

Table 2: Summary of LED ring animations indicating product malfunctions.

## 2.6 Battery life

Mode	Minimum working time	Conditions
Walking	11 h	Normal walking, level 4 damping

## 3. Mobile application

### 3.1 System requirements

The ENFORCE Prosthetic app is compatible with Android devices starting from version 10 (available on the Google Play store) and with iOS devices starting from version 15.5 (available on the App Store).

The application was tested on the following devices:

- Android: Motorola One Vision, Samsung Galaxy S20, Samsung Galaxy S22, Samsung Galaxy A54.
- iOS: iPhone 11 Pro, iPhone 11, iPhone X, iPhone Xs, iPhone8, iPhone 12 mini.

### 3.2 First connection to the prosthesis

**Attention:**

Before initiating the app for the first time, verify that your phone has Bluetooth enabled, and ensure the ENFORCE Bionic Foot is charged and functioning.

When launching the application for the first time, users are required to activate their user account. To proceed, click on the 'Activate account' button (1) on the login screen. Subsequently, enter your phone number (2) and the patient number (3), as specified in your contract with ENFORCE prosthetic center, then press the 'Activate account' button (4) once more.


Figure 6: Fill details.

In the final step, choose the preferred form of greeting displayed in the application (5), set a suitably strong password (6) and confirm it (7). Review the terms of use and privacy policy, agreeing to the conditions outlined (8). The activation process is concluded by pressing the 'Activate account' button (9).


Figure 7: Account activation.

### 3.3 Main screen

The home screen of the app enables users to monitor prosthesis performance and easily navigate essential functions.


Figure 8: Main screen.

1. User avatar - clicking on this opens a side menu.
2. Footwear Panel - displays the currently selected footwear.
3. Activity Panel - allows users to monitor daily step count and the percentage of their set goal achieved.
4. Battery Panel - enables users to check the current battery level.
5. Prosthesis Mode Panel - permits users to switch between different modes for the prosthesis.
6. Prosthesis Stiffness Panel - lets users adjust the stiffness of the damper.
7. Transition to My Prosthesis Screen.
8. Transition to Main Screen.
9. Transition to *Battery* Screen.

### 3.4 Side menu

Tapping on the user's avatar triggers the display of a side menu, simplifying navigation across different sections of the application.


Figure 9: Side menu.

1. My Account - users have the ability to modify personal data, profile photo, form of greeting, password, application language, and also delete their account.
2. My Prosthesis - users can connect to their device, update the prosthesis software, or add a measurement of the stump's circumference.
3. My Footwear - users can add, edit, delete, and calibrate footwear.
4. My Activity - users can monitor their activity statistics and adjust their daily step target as needed.
5. Notifications - users can access messages related to the operation of the prosthesis or application.
6. More - the user can read the regulations and policy privacy, find answers to frequently asked questions and submit your feedback or service request. Screen, more also allows you to use the Diagnostics functionality and the "Gestures and functions" screen.
7. Log Out - users can sign out of the application.

### 3.5 More screen

The screen allows you to use the prosthesis diagnostic functions, manage its gestures and functions. Additionally, you can use functions such as:


1. Service requests
2. Express your opinion about the application
3. Contact with a prosthetist
4. FAQ
5. Manual

Figure 10: More screen.

### 3.6 Connection to prosthesis

To initiate a connection with the prosthesis, navigate to the 'My Prosthesis' screen. This screen is accessible either from the side menu by clicking on the User Avatar (1) or by pressing the button with the Bluetooth symbol in the bottom menu (2), as illustrated in the *Figure 11a*.

a) Connecting with device.

b) Connecting with device step 2.

Figure 11: Connecting with prosthesis.

On the My Prosthesis screen, users can establish a connection with the prosthesis by clicking on its icon (1), as shown in *Figure 11b*. The color of the border indicates the current connection status between the prosthesis and the ENFORCE Prosthetic app:

- Red - not connected
- Blue - connected

Additionally, users can navigate to the Measurements screen (2), initiate an Update (3), or access the Details of their prosthesis (4) to view information such as the warranty period.

### 3.7 Automatic connection

After the device has been successfully paired with the mobile app for the initial setup, the app will automatically try to connect every time it is opened.

Please ensure that Bluetooth is turned on for both the mobile device and the electronic device to enable this automatic connection. You can toggle the automatic connection feature on or off by accessing the more section from the side menu.

### 3.8 Stump circumference measurements

In the Measurements screen, users can record stump circumference measurements once a day by clicking the Add button (1). The changes in circumferences are visualized in the 'Stump circumference statistics' graph (2). If you switch the graph display to a monthly date range, the results are presented as the average of the last 3 days. Additionally, the screen features a panel displaying measurements recorded on the current day (3). To edit the values, simply click the panel.


Figure 12: Stump circumference measurements.

### 3.9 Footwear calibration

On the My Footwear screen, the user can browse through the saved footwear. The currently selected footwear has a blue border, other footwear that has been added and calibrated has a grey border. Navigating between the saved items is done using the arrows (1 and 2).


Figure 13: Footwear management.

To add a new pair of footwear, press the Add Footwear button (3). Provide a name, and you can optionally add a photo or image from the phone's gallery. After clicking the Calibrate button, the calibration procedure will start, and it is essential to follow the on-screen instructions. The Delete Shoe button (4) allows you to remove an item, provided it is not the currently selected shoe. The Edit button (6) enables you to modify (change the image and name of) the currently displayed footwear. The Calibrate button (5) allows you to recalibrate. If the displayed footwear is changed, this button will be replaced by the Select button (5). Pressing it will adjust the settings of the prosthesis to match the selected footwear.

### 3.10 Reporting a problem

To access the Service Requests screen, tap on the Avatar, then choose More from the side menu, and finally, click on the Service requests button.

Figure 14: Service requests.

On the Service Requests screen, the user can check the details of their previous requests (1) or add new ones by clicking on the Add button (2).

### 3.11 Management of modes

The default operating mode of the prosthesis is walking mode. To go to the mode management panel, click the icon (5) in Figure 15. Then change the mode by clicking on the appropriate icon - bicycle, home or ladder mode (Figure 16). The currently selected option is highlighted in green. Changing to the default walking mode is possible by clicking again on the active mode (marked in green) - when no mode (bicycle, home or ladder) is highlighted in color, the prosthesis works in the default walking mode.


Figure 15: Main Screen.


Figure 16: Changing prosthesis mode.

### 3.12 Changing the stiffness of the prosthesis

To change the stiffness of the prosthesis, select the desired stiffness on the Main Screen (6, *Figure 17*), and follow the application prompts.


Figure 17: Main screen.

### 3.13 Features

This module allows you to turn on and off the functions and gestures of the prosthesis. By turning off a gesture or function, the user loses the ability to, for example, check the charge level of the prosthesis from the device.


Figure 18: Functions and gestures

1. Battery check. After turning it off, the gesture of checking the charge level of the prosthesis by tilting it to the maximum in plantar flexion will not be active. You will not be able to check the charge level on the LED Ring.
2. Bluetooth Management. After turning off the Bluetooth management function, the gesture to turn Bluetooth on or off will not be active, turning the prosthesis will not turn Bluetooth on or off, and the connection will be permanently turned on.
3. Home mode. After turning it off, the gesture of switching to home mode through maximum dorsiflexion is not active, which means that turning on home mode is only possible via the mobile application.
4. ARCA. Activation of active adaptation to the ground. The prosthesis still supports walking on varied terrain but adapts to it to a lesser extent.
5. Standing support - turning off this function will reduce the comfort when standing on the prosthesis and may extend the operating time on a single charge.

### 3.14 Software update

If a new software version becomes available after the application has been connected to the product, the user is presented with a dialog box informing them of the changes (*Figure 19*). After reviewing the changes, the user can initiate the update procedure by pressing the Update button (1) or postpone the update by closing the window with the "X" symbol (2).


Figure 19: Notification about available update and patch notes.

The mobile app and the device will then carry out the processes necessary for the update. The device will indicate that the update is ready by flashing the LED Ring pulsatingly, as described in *Table 1*. To proceed, press the Update button (1), as shown in *Figure 20*.


Figure 20: Begin update process.

When pressed, a system window will appear (*Figure 21*), prompting the user to confirm the connection to the Wi-Fi network created by the device. Click the Join button to proceed (1).

Figure 21: Connection to device's WiFi network.

The progress of the update transfer can be monitored in the window displayed by the application or directly on the LED Ring of the device, with the animation described in *Table 1*. Once the process is complete, the device will reset, and the app will display a message instructing you to reconnect to the device (*Figure 22*).


Figure 22: Complete update procedure.

### 3.15 Statistics and daily step goal

A daily step target can be set in the app. To do this, navigate to the statistics screen by selecting the item My Activity in the drop-down side menu (4, *Figure 9*) after clicking on the avatar (1, *Figure 8*). Another way to access this screen is by selecting a panel (3, *Figure 8*). While on the statistics screen, to edit a target, select the icon (1), as shown in *Figure 23a*. The screen in *Figure 23b* displays the current daily step target (2), the target change panel (3), and the save changes button (4).


Figure 23: Statistics screen (a) and steps target edit screen (b).

### 3.15 The self-diagnosis procedure

The self-diagnosis procedure allows you to check individual functionalities of the product. There are several tests to choose from:

1. The peripheral test allows you to determine whether all prosthesis components are functional.
2. Measurement test The measurement test allows you to determine whether the measurements made by the prosthesis sensors are correct.
3. The charging test allows you to determine whether the charging of the prosthesis is working properly, during this test you must connect the charger in accordance with the recommendations for charging the prosthesis.
4. The LED display test allows you to check whether all LEDs on the display are working properly. The display will change color during the test.
5. The damper test allows you to determine whether the prosthesis damper is working properly and has not been damaged.


Figure 24: The self-diagnosis procedure.

**Attention:**

Before initiating the app for the first time, verify that your phone has Bluetooth enabled, and ensure the ENFORCE Bionic Foot is charged and functioning.

After clicking the "Check" button, the testing procedure will begin. An information message will appear, after reading it, click the "Diagnose" button. This will start the test.

You will then receive a message about the success or failure of the test. If any of the tests fail, please contact us. Via "Service Requests". To submit a service request via the application, click on the Side Menu (*Figure 9*) and then on the "More" section (*Figure 10*).


Figure 25: The peripheral test.

#### 4. Conditions of use and storage

The product is IP68 waterproof which means it can be used in wet or humid environments and is resistant to contact with water. Wipe the product dry with a clean cloth after contact with fresh water or moisture.

If the product comes into contact with chlorinated water, chemicals, sand, dust or dirt, clean it immediately with fresh water and dry it. Ensure that the product is completely dry before reuse.

Temperature range for use of the product:  $-20^{\circ}\text{C}$  to  $+60^{\circ}\text{C}$ . Relative humidity up to 85%.

Storage and transport:

- 1 month temperature range of  $-20^{\circ}\text{C}$  to  $+60^{\circ}\text{C}$
- 3 months temperature range of  $-20^{\circ}\text{C}$  to  $+40^{\circ}\text{C}$
- 1 year temperature range of  $-20^{\circ}\text{C}$  to  $+20^{\circ}\text{C}$

If the product becomes dirty, proceed according to point 5. If the product is damaged, stop using it immediately and contact a prosthetics facility.

## 5. Maintenance and care

Keep the product clean. If dirt occurs, promptly wipe it off with a damp cloth. Pay particular attention to dirt and solids in the area of the magnetorheological damper (*Figure 1*, element 4), the connections of the upper and lower parts, and the charging socket.

**Attention:**

If the product shows signs of damage or wear, stop using it immediately.

## 6. Legal and safety notices

Signs of change or loss of utility function:

- Change in gait characteristics.
- Change in the position of components.
- Development of unusual noises in the components of the prosthesis.
- Leakage of fluid from the damper.

In order to maintain safety, it is necessary to:

- Use the product only during the authorized period of use set by the manufacturer (3 years).
- Undergo an annual safety check of the product.
- In the event of suspected damage to the product, check its function and serviceability.
- If the product is damaged, discontinue use immediately.
- Use only in permitted conditions and avoid harmful environments.
- The product must be checked for damage if it has been used in an unauthorized environment.
- Clean the contaminated product according to the instructions.
- When descending stairs and steep surfaces, use the handrail if available.
- It is strictly forbidden to use the product for extreme sports, running, bicycle racing, or other activities other than recreational cycling. These activities are undertaken at the user's own risk.
- Avoid contact of the charging cable with magnetic elements.
- Avoid contact of magnetic elements with the charging socket.

To maintain safety during the lifetime of the product, it is necessary to:

Description of Activities	Daily	Monthly	Annually
Check if the charging cable is damaged and if the charging socket is not contaminated	X		
Check if any fluid is leaking from the components of the prosthesis		X	
Check if the wires of the prosthesis are damaged		X	
Check the proper functioning of the prosthesis at the prosthetic facility			X
Check if the surface of the prosthesis is not greasy or contaminated with dust/sand	X		
Visual evaluation of the prosthesis		X	

## 7. Responsibility

The manufacturer (ENFORCE Medical Technologies LLC.) is liable in the event that the product is used in accordance with the recommendations contained in this document.

## 8. Warranty information

ENFORCE Medical Technologies provides a 24-month warranty on the product. Please refer to the product warranty card for details. Changes or modifications not expressly approved may void the warranty.

## 9. Electromagnetic compatibility

**Attention:**

Avoid using this appliance in the vicinity of or in combination with other equipment, as this may cause malfunction. If such use is necessary, inspect this unit and other equipment to verify proper operation.

**Attention:**

The use of accessories and cables other than those specified or supplied by the manufacturer of this equipment may result in increased electromagnetic emissions or reduced electromagnetic immunity of this equipment, leading to malfunctions.

**Attention:**

Portable RF communication equipment (including peripherals such as antenna cables and external antenna devices) should not be used closer than 30 cm (12 inches) to any part of the equipment, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment may occur. The device is suitable for use in any environment except flammable/explosive atmospheres or exposure to strong electric and/or magnetic fields (e.g., electrical transformers, high-power RF transmitters, RF surgical equipment, CT and MRI scanners).

### 9.1 Information on compliance regarding electromagnetic emissions

Test	Standard	Compliance level
Radiated RF emissions	<ul style="list-style-type: none"> <li>• Accordance with EN 55011:2016+A1:2017+A1:2020</li> <li>• Accordance with EN 606011- 2:2015-11/A1:2021</li> </ul>	Group 1- class B

### 9.2 Information on compliance regarding electromagnetic immunity

Test	Standard	Compliance level	
ESD resistance	EN 61000- 4- 2:2009	Indirect ESD:	+/- 8kV contact A
		Direct ESD:	+/- 8kV contact B +/- 15kV air B
Electromagnetic field immunity	EN 61000- 4- 3:2020	80 MHz - 2,7 GHz   10V/m   A 385MHz   27V/m   A 450 MHz   28V/m   A 710, 745, 780 MHz   9V/m   A 810, 870, 940 MHz   28V/m   A 1720, 1845, 1970, 2450 MHz   28V/m   A 5240, 5500, 5785 MHz   9V/m   A	
Immunity against radio fields in close distance	EN 61000-4-39:2017	8A/m 30 kHz /A 65 A/m 134.2 kHz /A 7,5 A/m 13.56 MHz /A	

## 10. Wireless modules

Technology	Frequency range	Power
Wi-fi	2412-2472 MHz / 2422-2462 MHz	19.19-19.53 dBm
Bluetooth LE	2402-2480 MHz	8.23Bm

## 11. Declaration of conformity (CE)

The product meets the requirements of Regulation (EU) 2017/745 of the European Parliament and of the Council of 5 April 2017 on medical devices.

The product complies with Directive 2014/30/EU of the European Parliament and of the Council of 26 February 2014 on the harmonisation of the laws of the Member States relating to electromagnetic compatibility.

## 12. Disposal

Dispose of the product and packaging in accordance with the current local or national environmental regulations.

Medical device

The device complies with Regulation (EU) 2017/745 of the European Parliament and of the Council, dated 5 April 2017, regarding medical devices.

ENFORCE Medical Technologies Sp. z o.o., ul. Szelałowska 45A 61-626 Poznań,  
Poland. Phone number: 61 111 00 69,  
email: [biuro@enforcemed.pl](mailto:biuro@enforcemed.pl), [www.enforcemed.pl](http://www.enforcemed.pl)

Serial number

The product is waterproof and dustproof

The product is equipped with a rechargeable battery, and it should not be disposed of with regular waste

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